

Product Datasheet

ARPC2 antibody (orb591894)

Description	Rabbit polyclonal antibody to ARPC2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of mouse ARPC2
Target	ARPC2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AVIH YR D D E T M Y V E S K K D R V T V V F S T V F K D D D D V V I G K V F M Q E F K E G R R A
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	Q9CVB6
NCBI	NP_083987.1
Expiration Date	12 months from date of receipt.