



Product Datasheet

Troponin I (phospho-Ser23/24) antibody (orb5920)


Description	Troponin I (phospho-Ser23/24) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Troponin I around the phosphorylation site of Ser23/24 (RR(p-S)(p-S)NY)
Target	Phospho-Cardiac Troponin I (Ser23 + Ser24)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	P19429
Dilution Range	Flow-Cyt=1ug/Test, WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human la...



Immunohistochemical staining of human ki...



Immunohistochemical staining of human la...