

Product Datasheet

PBLD2 antibody (orb592345)

Description	Rabbit polyclonal antibody to PBLD2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of mouse PBLD2
Target	PBLD2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EDPVTGSTHTLLGPYWSEELGKKEMRAFAQCSRRGGELDINLRPDGRVDIK
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	Q9CXN7
NCBI	NP_080361.1
Expiration Date	12 months from date of receipt.