


Product Datasheet

MSH2 antibody (orb592977)

Description	Rabbit polyclonal antibody to MSH2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human MSH2
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	150kDa
Uniprot ID	P43246
Entrez	4436
Dilution Range	WB: 1:500-1000, IF: 1:20-50, FACS: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeL...