


Product Datasheet

Cyclin E(phospho-Thr62) antibody (orb5930)

Description	Cyclin E(phospho-Thr62) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Cyclin E around the phosphorylation site of Thr77 (IP(p-T)PD)
Target	phospho-Cyclin E1 (Thr77)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	P24864
Dilution Range	Flow-Cyt=1µg /Test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of MCF7(blue line).Primar...