

Product Datasheet

TBR1 antibody (orb593013)

Description	Rabbit polyclonal antibody to TBR1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TBR1(NP_006584.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42kDa/74kDa
Uniprot ID	Q16650
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.