


Product Datasheet

[KO Validated] MYL6B antibody [Out of stock] (orb593102)

Description	Rabbit polyclonal antibody to [KO Validated] MYL6B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MYL6B (NP_002466.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	25kDa
Uniprot ID	P14649
Entrez	140465
Dilution Range	WB: 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from n...