

Product Datasheet

FBXW8 antibody (orb593119)

Description	Rabbit polyclonal antibody to FBXW8
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein of human FBXW8.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	67kDa
Uniprot ID	Q8N3Y1
Entrez	26259
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.