

Product Datasheet

SMG5 antibody (orb593120)

Description	Rabbit polyclonal antibody to SMG5
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human SMG5.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	114kDa
Uniprot ID	Q9UPR3
Entrez	23381
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.