

Product Datasheet

TBCCD1 antibody (orb593140)

Description	Rabbit polyclonal antibody to TBCCD1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human TBCCD1(NP_060608.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q9NVR7
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.