

Product Datasheet

SIX6 antibody (orb593158)

Description	Rabbit polyclonal antibody to SIX6
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human SIX6(NP_031400.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	095475
Dilution Range	IF 1:50-200
Expiration Date	12 months from date of receipt.