


Product Datasheet

PXN antibody (orb593167)

Description	Rabbit polyclonal antibody to PXN
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human Paxillin (NP_002850.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65KDa
Uniprot ID	P49023
Dilution Range	WB: 1:500-2000, IF: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...