

Product Datasheet

CFAP69 antibody (orb593178)

Description	Rabbit polyclonal antibody to CFAP69
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human CFAP69.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	A5D8W1
Entrez	79846
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.