

Product Datasheet

TBC1D4 antibody (orb593199)

Description	Rabbit polyclonal antibody to TBC1D4
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human TBC1D4(NP_055647.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160kDa
Uniprot ID	O60343
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.