

Product Datasheet

HSPB11 antibody (orb593205)

Description	Rabbit polyclonal antibody to HSPB11
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human HSPB11(NP_057210.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa
Uniprot ID	Q9Y547
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.