

Product Datasheet

ZHX2 antibody (orb593216)

Description	Rabbit polyclonal antibody to ZHX2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ZHX2(NP_055758.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120kDa
Uniprot ID	Q9Y6X8
Dilution Range	IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.