

Product Datasheet

Vgll4 antibody (orb593228)

Description	Rabbit polyclonal antibody to Vgll4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of mouse Vgll4(NP_808351.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31KDa
Uniprot ID	Q80V24
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.