

Product Datasheet

PPARA antibody (orb593232)

Description	Rabbit polyclonal antibody to PPARA
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human PPAR α .
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	52KDa
Uniprot ID	Q07869
Entrez	5465
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.