

Product Datasheet

CARTPT antibody (orb593253)

Description	Rabbit polyclonal antibody to CARTPT
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic peptide of human CARTPT(NP_004282.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12kDa
Uniprot ID	Q16568
Dilution Range	IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.