


Product Datasheet

DOCK2 antibody (orb593273)

Description	Rabbit polyclonal antibody to DOCK2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human DOCK2(NP_004937.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	240kDa
Uniprot ID	Q92608
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jurkat cells.