

Product Datasheet

NUP160 antibody (orb593285)

Description	Rabbit polyclonal antibody to NUP160
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human NUP160(NP_056046.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160kDa
Uniprot ID	Q12769
Dilution Range	WB: 1:500-1000, IHC-P: 1:20-50
Expiration Date	12 months from date of receipt.