


Product Datasheet

Tau (Phospho-Thr231) antibody (orb593322)

Description	Rabbit polyclonal antibody to Tau (Phospho-Thr231)
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T548 of human Tau(NP_058519.3?).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	60kDa
Uniprot ID	P10636
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse brain.