

Product Datasheet

MCL1 (Phospho-S159) antibody (orb593344)

Description Rabbit polyclonal antibody to MCL1 (Phospho-S159)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic phosphorylated peptide around S159 of human MCL1.

Preservatives Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Purity Affinity purification


MW 42kDa

Uniprot ID [Q07820](#)

Entrez [4170](#)

Dilution Range WB: 1:500-2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of 293...