

Product Datasheet

RAF1 (Phospho-S296) antibody (orb593350)

Description	Rabbit polyclonal antibody to RAF1 (Phospho-S296)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S296 of human Phospho-Raf1-S296.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P04049
Entrez	5894
Dilution Range	WB: 1:500-2000, IP: 1:50-100
Expiration Date	12 months from date of receipt.