

Product Datasheet

Inactivated Varicella Zoster Virus protein (orb593477)

Description	Inactivated Varicella Zoster Virus protein
Species/Host	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized powder or 0.1M PBS, PH 7.0
Storage	Store the lyophilized powder at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Purity	> 95%
Source	Virus
Expiration Date	6 months from date of receipt.