

Product Datasheet

Cytomegalovirus antibody (orb593536)

Description	Mouse monoclonal antibody to Cytomegalovirus
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Inactivated Cytomegalovirus
Preservatives	0.1M PBS, PH 8.0
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.