

## Product Datasheet

# Pneumococcal Type 12F Polysaccharide antibody (orb593567)

|                                |                                                                                                                                                                              |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Rabbit polyclonal antibody to<br>Pneumococcal Type 12F<br>Polysaccharide                                                                                                     |
| <b>Reactivity</b>              | Rabbit                                                                                                                                                                       |
| <b>Conjugation</b>             | Unconjugated                                                                                                                                                                 |
| <b>Tested<br/>Applications</b> | ELISA, Nephelometry, WB                                                                                                                                                      |
| <b>Immunogen</b>               | Pneumococcal Type 12A<br>Polysaccharide                                                                                                                                      |
| <b>Preservatives</b>           | 0.1M PBS, PH 8.0                                                                                                                                                             |
| <b>Storage</b>                 | Store at 4°C for 6 months. It is<br>recommended to aliquote to ~10ul<br>per vial and store at -20°C for long<br>term storage. Avoid repeated<br>freezing and thawing cycles. |
| <b>Note</b>                    | For research use only                                                                                                                                                        |
| <b>Expiration Date</b>         | 12 months from date of receipt.                                                                                                                                              |