

## Product Datasheet

### Pneumococcal Type II Polysaccharide antibody (orb593578)

|                            |                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to<br>Pneumococcal Type II Polysaccharide                                                                                             |
| <b>Reactivity</b>          | Rabbit                                                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                     |
| <b>Tested Applications</b> | ELISA, Nephelometry, WB                                                                                                                                          |
| <b>Immunogen</b>           | Pneumococcal Type II Polysaccharide                                                                                                                              |
| <b>Preservatives</b>       | 0.1M PBS, PH 8.0                                                                                                                                                 |
| <b>Storage</b>             | Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles. |
| <b>Note</b>                | For research use only                                                                                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                  |