

Product Datasheet

Pneumococcal Type 23F Polysaccharide antibody (CRM197) (orb593581)

Description	Rabbit polyclonal antibody to Pneumococcal Type 23F Polysaccharide (CRM197)
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, Nephelometry, WB
Immunogen	Pneumococcal Type 23F Polysaccharide CRM197 Conjugate
Preservatives	0.1M PBS, PH 8.0
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.