

Product Datasheet

Pneumococcal Type VIII Polysaccharide antibody (orb593595)

Description	Rabbit polyclonal antibody to Pneumococcal Type VIII Polysaccharide
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, Nephelometry, WB
Immunogen	Pneumococcal Type VIII Polysaccharide
Preservatives	0.1M PBS, PH 8.0
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.