

Product Datasheet

Sabin Poliovirus Type III antibody (orb593603)

Description	Rabbit polyclonal antibody to Sabin Poliovirus Type III
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Inactivated Sabin Poliovirus Type III
Preservatives	0.1M PBS, PH 8.0
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.