

Product Datasheet

Diphtheria Toxoid antibody (HRP) (orb593633)

Description	Goat polyclonal antibody against Diphtheria Toxoid conjugated to HRP
Reactivity	Goat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified Diphtheria Toxoid
Preservatives	0.1M PBS, PH 8.0, with protein stabilizer
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.