

Product Datasheet

Vero Cell Protein (HCP) antibody (HRP) (orb593650)

Description	Rabbit polyclonal antibody against Vero Cell Protein (HCP) conjugated to HRP
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Vero Cell Lysate
Preservatives	0.1M PBS, PH 8.0, with protein stabilizer
Storage	Store at 4°C for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.