

Product Datasheet

CYP11B1 antibody (orb5937)

Description

CYP11B1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC-P, WB

Immunogen

KLH conjugated synthetic peptide derived from human CYP11B1. Please contact us for the exact immunogen sequence. The peptide is available as orb12918.

Target

CYP11B1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

58 kDa

Uniprot ID

[P15538](#)

NCBI

[005250864](#), [005250807](#), [11](#)

Entrez

[1584](#)

Dilution Range


IF/ICC: 1:50-500, IHC-P: 1:50-300, WB: 1:200-2000

Expiration Date


12 months from date of receipt.



WB analysis of mouse brain (lane 1), rat...



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of paraffin...