


Product Datasheet

Mouse Csf1 protein (orb594718)

Description	Recombinant Mouse Macrophage colony-stimulating factor 1(Csf1),partial (Active)
Endotoxins	Less than 0.01 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	KEVSEHCSHMIGNGHLKVLQQLIDSQMETSCQIAFEFVDQEQQLDDPVCYLKKAF FLVQDIIDETMRFKDNTNPNANATERLQELSNLNSCFTKDYEEQNKACVRTFHET PLQLEKIKNFFNETKNLLEKDWNIFTKNCNNSFAKCSSRDVVTKPDCNCLYPKA TPSSDPASASPHQPPAPSMAPLAGLAWDDSQRTGESSLLPSELPLRIEDPGSAKQ RPPRSTCQTLE
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	26 kDa
Uniprot ID	P07141
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous
SDS-PAGE...