

Product Datasheet

Human PDGF B protein (orb594736)

Description	Recombinant Human Platelet-derived growth factor subunit B(PDGFB) (Active)
Endotoxins	Less than 0.2 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered solution of 20mM NaAc-HAc, pH 4.5.
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	SLGSLTIAEPAMIAECKTRTEVF FEIS RRRLIDRTNANFLVWPPCVEVQRCSGCCNNR NVQCRPTQVQLRPVQVRKIEIVRKKPIFKKATVTLEDHLACKCETVAAARPVT
Purity	Greater than 98% as determined by SDS-PAGE.
Source	E.coli
MW	12.42 kDa
Uniprot ID	P01127
Expiration Date	6 months from date of receipt.