


Product Datasheet

Human IL22 protein (orb594790)

Description	Recombinant Human Interleukin-22(IL22),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 20mM Histidine-HCl, 6% Sucrose, 4% Mannitol, 0.05% Tween 80, pH 5.5.
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	APISSHCRLDKSNFQQPYITNRTFMLAKEASLADNNTDVRLLIGEKLFGVSMSE CYLMKQVLNFTLEEVLPQSDRFQPYMQEVVPFLARLSNRLSTCHIEGDDLHIQR NVQKLKDTVKKLGESGEIKAIGELDLLFMSLRNACI
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	16.9 kDa
Uniprot ID	Q9GZX6
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous
SDS-PAG...