

Product Datasheet

Human IL23R protein (orb594791)

Description	Recombinant Human Interleukin-23 receptor(IL23R),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal Fc-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	GITNINCSGHIWVEPATIFKMGMNISIYCQAAIKNCQPRKLHFKNGIKERFQITRI NKTARLWYKNFLEPHASMYCTAECPKHFQETLICGKDISSGYPPDIPDEVTCVIY EYSGNMTCTWNAGKLTYIDTKYVVHVKSLETEEQQYLTSSYINISTDSLQGGKK YLVVWQAANALGMEESKQLQIHLDDIVIPSAAVISRAETINATVPKTIYWDSQTTI EKVSCEMRYKATTNQTWNVKEFDTNFTYVQQSEFYLEPNIKYVFQVRCQETGKR YWQPWSSLFFHKTPETVPQVTSKAFQHDTWNSGLTVASISTGHLTSDNRGD
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	65 kDa
Uniprot ID	Q5VWK5
Expiration Date	6 months from date of receipt.