


Product Datasheet

Mouse IL4 protein (orb594835)

Description	Recombinant Mouse Interleukin-4(IL4),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 300mM NaCl, 5% Trehalose, pH 6.5.
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	HGCDKNHLREIIGILNEVTGEGTPCTEMDVPNVLTATKNTTESELVCRASKVLRIFYLKHGKTPCLKKNSSVLMELQRLFRAFRCCLDSSISCTMNESKSTSLKDFLESLSIMQMDYS
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	13.4 kDa
Uniprot ID	P07750
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous
SDS-PAG...