


Product Datasheet

Human RANKL protein (orb594839)

Description	Recombinant Human Tumor necrosis factor ligand superfamily member 11(TNFSF11),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 20 mM Tris-HCl, 150 mM NaCl, pH 8.0
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag free
Note	For research use only
Application notes	Partial
Protein Sequence	IRAEKAMVDGSWLDLAKRSKLEAQPF AHLTINATDIPSGSHKVSLSWYHDRGW AKISNMTFSNGKLIVNQDGFYYLYANICFRHHETSGDLATEYLQLMVYVTKTSIKIP SSHTLMKGGSTKYWSGNSEFHFYSINVGGFFKLRSGEISIEVSNPSLLDPDQDA TYFGAFKVRDID
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	22.4 kDa
Uniprot ID	O14788
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous
SDS-PAGE...