


Product Datasheet

Human TNF protein (orb594840)

Description	Recombinant Human Tumor necrosis factor(TNF),partial (Active)
Endotoxins	Less than 0.01 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 6% Sucrose, 4% Mannitol, 0.05% Tween 80, pH 6.0
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	VRSSSRTPSPDKPVAHVVANPQAEGQLQWLNRRANALLANGVELRDNQLVVPSE GLYLIYSQVLFGQGCPSTHVLLTHTISRIAVSYQTKVNLLSAIKSPCQRETPEGAE AKPWYEPIYLGGVFQLEKGDRLSAEINRPDYLDFAESGQVYFGIIAL
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	17.5 kDa
Uniprot ID	P01375
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous
SDS-PAGE...