

Product Datasheet

Human TNFSF10 protein (orb594844)

Description	Recombinant Human Tumor necrosis factor ligand superfamily member 10(TNFSF10),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	0.2 μm filtered 40mM Tris-HCl, 300mM NaCl, 5%Trehalose, 5%Mannitol, 0.01%Tween80, 10%Glycerol, pH 7.0.
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Partial
Protein Sequence	RERGPQRVAAHITGTRGRSNTLSSPNSKNEKALGRKINSWESSRSGHSFSLNLH LRNGELVIHEKGFYIYSQTYFRFQEEIKENTKNDKQMVQYIYKYTSYPDPILLMKS ARNSCWSKDAEYGLYSIQGGIFELKENDRIFVSVTNEHLIDMDHEASFFGAFLV G
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	19.5 kDa
Uniprot ID	P50591
Expiration Date	6 months from date of receipt.

