

Product Datasheet

Human TNFRSF10B protein (orb594847)

Description	Recombinant Human Tumor necrosis factor receptor superfamily member 10B(TNFRSF10B),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 20 mM PB, 150 mM NaCl, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	ITQQDLAPQQRAAPQQRSSPSEGLCPPGHHISEDGRDCISCKYGQDYSTHWN DLLFCLRCTRCDSGEVELSPCTTTRNTVCQCEEGTFREEDSPEMCRKCRTGCPR GMVKVGDCTPWSDIECVHKE
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	15.19 kDa
Uniprot ID	O14763
Expiration Date	6 months from date of receipt.