

Product Datasheet

Human TNFRSF9 protein (orb594867)

Description	Recombinant Human Tumor necrosis factor receptor superfamily member 9(TNFRSF9),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-Fc-tagged
Note	For research use only
Application notes	Extracellular Domain
Protein Sequence	LQDPCSNCPAGTFCDNNRNQICSPCPPNSFSSAGGQRTCDICRQCKGVFRTTRKE CSSTSNAECDCTPGFHCLGAGCSMCEQDCKQGQELTKKGCKDCCFGTFNDQK RGICRPWTNCSLDGKSVLVNGTKERDVVCGPSPADLSPGASSVTPPAPAREPGH SPQ
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	44 kDa
Uniprot ID	Q07011
Expiration Date	6 months from date of receipt.