

Product Datasheet

Human CA8 protein (orb594896)

Description	Recombinant Human Carbonic anhydrase-related protein(CA8),partial (Active)
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	0.2 μm filtered 20 mM Tris-HCl, 500 mM NaCl, 1 mM DTT, pH 8.5
Form/Appearance	Liquid
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	ADLSFIEDTVAFPEKEEDEEEEEEGVEWGYEEGVEWGLVFPDANGEYQSPINLN SREARYDPSLLDVRLSPNYVVCRDCEVTNDGHTIQVILKSKSVLSGGPLPQGHEF ELYEVRFHWGRENQRGSEHTVNFKAFPMELHLIHWNSTLFGSIDEAVGKPHGIAI IALFVQIGKEHVGLKAVTEILQDIQYKGKSKTIPCFNPNTLLPDLLRDYWVYEGSL TIPPCSEGVTWILFRYPLTISQLQIEEFRRRLRTHVKGAELVEGCDGILGDNFRPTQP LSDRVIRAAFAQ
Purity	Greater than 95% as determined by SDS-PAGE.
Source	E.coli
MW	34.04 kDa
Uniprot ID	P35219
Expiration Date	6 months from date of receipt.