

Product Datasheet

Human CILP protein (orb594963)

Description	Recombinant Human Cartilage intermediate layer protein 1(CILP),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal MBP-tagged and C-terminal 6xHis-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	EDRTFLVGNLEIRERRLFNLDVPESRRRCFVKVRAYRSERFLPSEQIGGVVISVINLE PRTGFLSNPRAWGRFDSVITGPNACVPAFCDDQSPDAYSAYVLASLAGEELQA VESSPKFNPNAIGVPQPYLNKLNRYRTDHEDPRVKKTAFQISMAKPRPNSAEESN GPIYAFENLRACEEAPPSAAHFRFYQIEGDRYDYNTPFNEDDPMSWTEDYLAW WPKPMEFRACYIKVKIVGPLEVNVRSRNMGGTHRQTVGKLYGIRDVRSTRDRD QPNVSAACLEFKCSGMLYDQDRVDRTLKVIPQGSCRRASVNPMLHEYLVNHLP LAVNNDTSEYTMLAPLDPLGHNYGIYTVTDQDPRTAKEIALGRCFDGTSDGSSRI MKSINVGVALTFNCVERQVGRQSAFQYLQSTPAQSPAAGTVQGRVPSRRQRA SRGGQRQGGVVASLRFPVRAQQPLIN
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	95.7 kDa
Uniprot ID	075339
Expiration Date	6 months from date of receipt.