

Product Datasheet

Human PSMD14 protein (orb595011)

Description	Recombinant Human 26S proteasome non-ATPase regulatory subunit 14(PSMD14)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal MBP-tagged and C-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length
Protein Sequence	MDRLLRLGGGMPGLGQGPPTDAPAVDTAEQVYISSLALLKMLKHGRAGVPMEV MGLMLGEFVDDYTVRVIDVFAMPQSGTGVSV EAVDPVFQAKMLDMLKQTRPE MVVGWYHSHPGFGCWLSGVDINTQQSFEALSERAVAVVVDPIQSVKGKVVIDA FRLINANMMVLGHEPRQTTSNLGHLNKPISQALIHGLNRHYYSITINYRKNELEQK MLLNLHKKSWMEGLTLQDYSEHCKHNESVVKEMLELAKNYNKAVEEEDKMTPE QLAIKNVGKQDPKRHLEEHVDVLMTSNIVQCLAAMLDTVVFK
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	79.2 kDa
Uniprot ID	O00487
Expiration Date	6 months from date of receipt.