

Product Datasheet

Primate TMEM168 protein (orb595030)

Description	Recombinant Pongo abelii Transmembrane protein 168(TMEM168),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	EWWRKNGSFC SRLIIVLDSENSTPWVKEVRKINDQYIAVQGAELIKTVDIEEADP PQLGDFTKDWVEYNCNSSNNICWTEKGRTVKAVYGVSKRWSDYTLHLPTGSDV AK
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	16.8 kDa
Uniprot ID	Q5RD28
Expiration Date	6 months from date of receipt.