

Product Datasheet

Mouse Fam20c protein (orb595090)

Description	Recombinant Mouse Extracellular serine/threonine protein kinase FAM20C(Fam20c)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	DFSSDPASNLTSHSLEKLPSAAEPVDHAPRGQEPRSPPPRDPAPHRPLLRDPGPRP RVPPPGPSGDGSLAKLFEHPLYQGAVPPLTEDDVLFNVNSDIRFNPKAAENPD WPHEGAEGAEFLPTGEAAVNLYPNWLKFHIGINRYELYSRHNPDAIDALLRDLGSQ KITSVAMKSGGTQLKLIMTFQNYGQALFKPMKQTREQETPPDFFYFSDYERHNAE IAAFHLDRILD FRRVPPVAGRMINMTKEIRDVTRDKKLWRTFFVSPANNICFYGEC SYYCSTEHALCGRPDQIEGSLAAFLPDLSLAKRKTWRNPWRRSYHKRKAWEV DPDYCEEVKQTPPYDSGHRILDIMDMTVFDFLMGNMDRHHYETFEKFGNETFIIH LDNGRGFGKYSHDELSILAPLHQCCRIRRSTYLRLLAKEEHKLSLLMAESLQH DKVAPVLYQLHLEALDRRLRIVLQAVRDCVEKDGLSSVVEDDLATEHRASTER
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	60.4 kDa
Uniprot ID	Q5MJS3
Expiration Date	6 months from date of receipt.