



Product Datasheet

CYP2C9 antibody (orb5951)

Description	CYP2C9 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CYP2C9 (351-450/490 aa)
Target	CYP2C9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P11712
Dilution Range	WB=1:500-2000, IHC-P=1:200-800
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of rat liver...